



Flame
Red
Frame



Flame
Blue
Frame

NEW!




**EASY HEIGHT
ADJUSTMENT
WITH SELF
THREADING
KNOB!**



10289RD



**COMES WITH CARRY
BASKET AND TRAY!**



815B

**BASKET CONVENIENTLY
FOLDS FOR EASY STORAGE!**

Deluxe 3 Wheel Aluminum Rollator, 8" Casters with Loop Locks

Available
1 Pack

10289RD Flame Red, (Has replaced 10218), 1/cs
10289BL Flame Blue, (Has replaced 10218), 1/cs
815B Basket, 1/ea

FEATURES:

- Comes standard with basket, tray and pouch.
- Ships in ultra compact carton. Comes with frame folded and compact. Simply loosen lock knob and pull trombone style frame to lock in place.
- Lightweight aluminum frame.
- Lightweight, solid 8" wheels for indoor or outdoor use.
- Adjustable handle height.
- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- Easy, one hand folding.
- Limited Lifetime Warranty.

For most updated parts diagrams, please visit our website
@ www.drivemedical.com

Winnie Lite Supreme/Go Lite 3 Wheel Aluminum Rollators, 8" Casters and Loop Locks

Available
1 Pack

199 Tan Plaid, 1/cs
10218FBL-1* Flame Blue, (Has replaced 196), 1/cs
10218FRD-1* Flame Red, (Has replaced 195), 1/cs

FEATURES:

- Special loop lock made of internal aluminum casting operates easily and ensures safety.
- 8" casters with soft-grip tires are ideal for indoor and outdoor use.
- Comes standard with carry pouch.
- Handles are adjustable in height.
- Easy, one hand folding.
- Optional wire basket available (item # 806).
- Optional basket ring available (item #160).
- Limited Lifetime Warranty.




10218FRD-1



**EASY HEIGHT
ADJUSTMENT
WITH SELF
THREADING
KNOB!**



Tan Plaid
Bag



Tan
Frame



Flame
Red
Frame



Flame
Blue
Frame



199



**PARTS
DIAGRAMS
ON PAGE 71 & 72**

SPECIFICATIONS

Item #	Handle Height (highest)	Handle Height (lowest)	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
10289 SERIES	38"	31"	25"	24"	NA	NA	NA	16.5 lbs.	20.5 lbs.	300 lbs.
WINNIE LITE SUPREME	38"	32"	26.5"	24"	NA	NA	NA	12 lbs.	16 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

rollators

drive

To order call toll free: 877.224.0946