

SMALL BASE QUAD CANES CONTINUED





**THIS ITEM
IS AVAILABLE
IN A 1-PACK!**

Quad Canes, Small Base, Designer Colors

10312CK-1	Foam Rubber Grip, Checker Cane, (Has replaced 10312CK-4), 1/cs
10312FP-1	Foam Rubber Grip, Floral Design Cane, (Has replaced 10312FB-4), 1/cs
10312RD-1	Foam Rubber Grip, Rainbow Cane, (Has replaced 10312RD-4), 1/cs
10312SF-1	Foam Rubber Grip, Nautical Cane, (Has replaced 10312SF-4), 1/cs
10312BE-4*	Foam Rubber Grip, Birds Eye Cane, 4/cs
10312FP-4*	Foam Rubber Grip, Floral Design Cane, 4/cs

FEATURES:

- 4 point base, combined with offset handle provides additional stability and support.
- Easy-to-use, one-button height adjustment with locking bolt, allows quick, safe fitting.
- Height adjustment 30" - 39".
- No assembly required.
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.