



THIS ITEM
IS AVAILABLE
IN A 1-PACK!



NEW
TAB
SILENCER



PREVENTS
RATTLING!

10381SL-6
Silver


10381BLK-6
Black



10349



10349HCTR



MEDICARE
REIMBURSEMENT CODE
E0105
WHEN ASSEMBLED
TO ANY CANE

Offset Aluminum Cane with Tab-Loc Silencer

10381SL-6	Silver, 6/cs
10381SL-1	Silver, 1/cs
10381BLK-6	Black, 6/cs
10381BLK-1	Black, 1/cs
10381BLKXL-6	Black, Extra Long, 6/cs

FEATURES:

- Tab-Loc silencer provides added safety to ensure a locked height adjustment.
- Features a padded grip. The flat side of the grip faces up to reduce stress that the user's hand may experience.
- Ergonomically designed padded handle.
- Comes standard with wrist strap.
- Manufactured with an aluminum shaft and a steel base to provide the perfect balance of safety and durability.
- Handle height adjusts from 29" to 38".
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.

Tri-Support Cane Tip

10349HCTR For Use with 3/4" Tubing, 1/bx

FEATURES:

- Replaces existing cane tip.
- Provides better balance and stability than a conventional cane tip.
- Allows cane to stand on its own.
- Three tier design allows for better traction.
- For use with 3/4" tubing.
- Fits most manufacturers canes.

SEE
PAGE 107
FOR SIZING
CHART

AbleTripod® Cane Tip

10349 Fits 3/4" Tubing, 1/ea

FEATURES:

- Replaces existing cane tip.
- Provides balance, stability and makes your cane self-standing.
- For right or left handed base.
- Flexes as you walk reducing impact through hand & wrist.
- Easy to use on stairs and provides greater traction on most terrain including grass, gravel and sand.
- Fits most manufacturers canes.

CANE
DISPLAY
AVAILABLE

SEE PAGE 102

CAN BE USED WITH ALL DRIVE OFFSET, PALM GRIP, "T" HANDLE, FOLDING, CROOK. ROUND HANDLE AND WOODEN CANES!

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.