



**THIS ITEM
IS AVAILABLE
IN A 1-PACK!**

Palm Grip Cane (Fisher Style™)

10339R-6	Right, 6/cs
10339L-6	Left, 6/cs
10339R-1	Right, 1/cs
10339L-1	Left, 1/cs

FEATURES:

- Unique "Palm Grip" is designed for right handed person or left handed person.
(Specify when ordering)
- Easy-to-use, one-button height adjustment with locking ring ensures added safety after adjustment.
- Height adjusts from 29" - 38".
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.



**THIS ITEM
IS AVAILABLE
IN A 1-PACK!**

Blind Folding Cane

10352-1 (Has replaced 10352-6), 1/cs

FEATURES:

- High quality 4 section aluminum construction.
- Shaft covered with white and red reflective tape for night visibility.
- Cane length: 46 1/4".
- Comes standard with strap and reinforced nylon tip.

