



**THE SLING SEAT
EASILY FOLDS FOR
CONVENIENT
TRANSPORT AND OPENS
TO REST ON WHEN
NEEDED!**



10360

Cane/Sling Seat

10360

1/cs

FEATURES:

- Provides a comfortable seat to rest on when open, and a sturdy quad cane when closed.
- Strong aluminum bronze construction.
- Four sturdy legs with vinyl contoured tips.
- Nylon sling bench, 18" x 8".
- Foam grip handles.
- 34" tall when closed.
- **Seat to floor height: 22"**
- Weighs 2.4 lbs.
- **250 lb. Weight Capacity.**
- Limited Lifetime Warranty on frame, 6 Month Warranty on Sling.



**NEW
RETAIL
PACKAGE!**



Cane Seat

10365

Seat to Floor Height: 20", 1/cs

10365-ADJ

Adjustable Height

Cane Height: 34"-38",

Seat Height: 19"-22", 1/cs



**NEW
RETAIL
PACKAGE!**



FEATURES:

- Provides a comfortable seat to rest on when open, and a sturdy support cane when closed.
- Strong aluminum construction.
- Tri-pod design with vinyl contoured tipped legs.
- 8.4" diameter seat, 34" tall when closed.
- **Seat to floor height: 20"**
- Weighs 1.3 lbs.
- **250 lb. Weight Capacity.**
- Limited Lifetime Warranty.



10365

**THE CANE SEAT
EASILY FOLDS FOR
CONVENIENT
TRANSPORT AND
OPENS
TO REST ON
WHEN NEEDED!**



10365-ADJ

canes & crutches

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.