




THIS ITEM
IS AVAILABLE
IN A 1-PACK!


COMPARE TO:
drive Invacare®
10385-4 7827



NEW TAB
SILENCER
PREVENTS
RATTLING!

NEW SILVER VEIN FINISH!
LESS INSTITUTIONAL &
EASY TO MAINTAIN!

ADDITIONAL
QUAD CANES
ON PAGE 72-79



10385-1



Quad Cane, Large Base, with Silver Vein Finish

10385-1 (Has replaced 10300-1), 1/cs
10385-4 (Has replaced 10385-2), 4/cs


FEATURES:

- Tab-Loc silencer provides added safety to ensure a locked height adjustment.
- Features a padded grip. The flat side of the grip faces up to reduce stress that the user's hand may experience.
- 4 point base, combined with offset handle provides additional stability and support.
- Attractive silver vein finish on base.
- Ergonomically designed padded handle.
- Manufactured with an aluminum shaft and a steel base to provide the perfect balance of safety and durability.
- Handle height adjusts from 29" to 38".
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.



THIS ITEM
IS AVAILABLE
IN A 1-PACK!


COMPARE TO:
drive Invacare®
10384-4 3920-2



NEW TAB
SILENCER
PREVENTS
RATTLING!

NEW SILVER VEIN FINISH!
LESS INSTITUTIONAL &
EASY TO MAINTAIN!

ADDITIONAL
QUAD CANES
ON PAGE 72-79



10384-1



Quad Cane, Small Base with, Silver Vein Finish

10384-1 (Has replaced 10301-1), 1/cs
10384-4 (Has replaced 10384-2), 4/cs

FEATURES:

- Tab-Loc silencer provides added safety to ensure a locked height adjustment.
- Features a padded grip. The flat side of the grip faces up to reduce stress that the user's hand may experience.
- 4 point base, combined with offset handle provides additional stability and support.
- Attractive silver vein finish on base.
- Ergonomically designed padded handle.
- Manufactured with an aluminum shaft and a steel base to provide the perfect balance of safety and durability.
- Handle height adjusts from 29" to 38".
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.



THIS ITEM
IS AVAILABLE
IN A 1-PACK!

HANDLE
EASILY
ADJUSTS
FOR LEFT
OR RIGHT
HAND USE!



REINFORCED
STEEL
PLATE

10317-4



10316-4

Bariatric Quad Canes

10316-4 Small Base, Vinyl Contoured Hand Grip, Chrome Cane, 4/cs
10317-1 Large Base, Vinyl Contoured Hand Grip, Chrome Cane, 1/cs
10317-4 Large Base, Vinyl Contoured Hand Grip, Chrome Cane, 4/cs

FEATURES:

- 500 lb. Weight Capacity.
- Lightweight, extra strong base.
- Easy-to-use, one button height adjustment from 30" - 39".
- Manufactured with 1" aluminum tubing, steel base.
- Comfortable vinyl, contoured hand grips.
- Aluminum handle, steel base.
- Limited Lifetime Warranty.

BARIATRIC
500 LB
WEIGHT CAPACITY

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

canes & crutches

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM