

10410
Silver10411*
Red10412
Blue10413*
Green10414*
Pink10407B*
Blue10407R*
Red10403G*
Chrome10403
Chrome

BARIATRIC
500 LB
WEIGHT CAPACITY

EuroStyle Lightweight Aluminum Forearm Crutch

10410	Adult, Silver, 1 pr/bx
10411*	Adult, Red, 1 pr/bx
10412	Adult, Blue, 1 pr/bx
10413*	Adult, Green, 1pr/bx
10414*	Adult, Pink, 1 pr/bx

FEATURES:

- One piece molded plastic cuff and hand grip assembly provides safety and comfort.
- Available in attractive colors.
- Limited Lifetime Warranty.

Steel, Forearm Crutches

10403	Adult,
10404*	Youth, 1 pr/bx
10405	300 lb. Weight Capacity, 1pr/bx Tall Adult,
10407B*	300 lb. Weight Capacity, 1 pr/bx Child, Blue, 1 pr/bx
10407R*	Child, Red, 1 pr/bx
10403G*	Ortho K Grip, Adult, 1 pr/bx
10404G*	Ortho K Grip, Youth, 1 pr/bx
10405G*	Ortho K Grip, Tall Adult, 1 pr/bx

Bariatric Steel, Forearm Crutches

10403 HD	Adult, Bariatric,
10404 HD*	500 lb. Weight Capacity, 1pr/cs Youth, Bariatric,
10405 HD*	500 lb. Weight Capacity, 1pr/cs Tall Adult, Bariatric,
	500 lb. Weight Capacity, 1pr/cs

FEATURES:

- Vinyl coated, ergonomically contoured arm cuffs molded for comfort and stability.
- Leg and forearm sections adjust independently for optimal sizing.
- Vinyl hand grips are comfortable and durable.
- Extra large tips for added stability.
- Ortho K handle is hard plastic.
- Limited Lifetime Warranty.

SPECIFICATIONS

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Item #	Handle to Floor (highest setting)	Handle to Floor (lowest setting)	Patient Height (highest)	Patient Height (lowest)	Product Weight (each)	Weight Limit
10410 - 10414	39"	30"	6'2"	5'0"	1.5 lbs.	300 lbs.
10403/G	37"	28"	6'2"	5'0"	2.5 lbs.	300 lbs.
10404/G	30.5"	21.5"	5'2"	4'2"	2.3 lbs.	300 lbs.
10405/G	42"	33"	6'6"	5'10"	2.6 lbs.	300 lbs.
10407B/R	21"	14"	4'6"	3'2"	.93 lbs.	250 lbs.
10403 HD	37"	28"	6'2"	5'0"	5.7 lbs.	500 lbs.
10404 HD	30.5"	21.5"	5'2"	4'2"	4.6 lbs.	500 lbs.
10405 HD	42"	33"	6'7"	5'10"	6.8 lbs.	500 lbs.