



10407R



10404G



HPCPS
E0110 / E0111

Aluminum Forearm Crutches

Item #	Product Handgrip	Product Size	Product Handle to Floor Height	Patient Height	Product Weight	Product Carton Shipping Weight	Product Weight Limit	UOM
10407R	Vinyl	Child	14" - 21"	3'2" - 4'2"	1 lb.	3 lbs.	250 lbs.	1 pr/bx
10404G	Hard Plastic	Youth	21.5" - 30.5"	4'6" - 5'2"	2.3 lbs.	5 lbs.	300 lbs.	1 pr/bx

- Vinyl-coated, ergonomically contoured arm cuffs molded for comfort and stability.
- Leg and forearm sections adjust for optimal sizing.
- Extra-large tips for added stability.