



10400-1



10425-8

## Aluminum Crutches with Accessories

10400-1	Adult, 1 pr/cs
10400-8	Adult, 8 pr/cs
10401-1	Youth, 1 pr/cs
10401-8	Youth, 8 pr/cs
10402-1	Tall Adult, 1 pr/cs
10402-8	Tall Adult, 8 pr/cs
10416-1	Pediatric, (Has replaced 10416-2), 1 pr/cs



**THIS ITEM  
IS AVAILABLE  
IN A 1-PACK!**

10400C-8*	Combo Pack, 8pr/cs
<b>COMBO PACK</b>	Contains 4 pr Adult, 2 pr Tall Adult Adult, 2 pr Youth

### FEATURES:

- Lighter than wood, the standard aluminum crutch provides stability and durability.
- Double extruded center tube provides additional strength to weight-bearing area.
- Push-pin adjustment makes it easy to adjust crutch height in 1" increments.
- Non-skid, jumbo size vinyl contoured tips provide excellent traction.
- Easy wing nut adjustments for hand grips.
- Comfortable, durable underarm pad and hand grip.
- 350 lb. Weight Capacity on aluminum crutch (excluding pediatric).
- Limited Lifetime Warranty.

## Wooden Crutches with Accessories

10425-1	Adult, 1 pr/cs
10425-8	Adult, 8 pr/cs
10426-1	Youth, 1 pr/cs
10426-8	Youth, 8 pr/cs
10427-1	Tall Adult, 1 pr/cs
10427-8	Tall Adult, 8 pr/cs



**THIS ITEM  
IS AVAILABLE  
IN A 1-PACK!**

### FEATURES:

- Laminated wood crutches are made from the finest hardwood veneers.
- Crutches come complete with accessories (closed hand grips, underarm cushions and tips).
- Easy wing-nut adjustment for height and hand grips.
- 250 lb. Weight Capacity.
- Limited Lifetime Warranty.

### SPECIFICATIONS

Item #	Height (highest setting) (armrest to floor)	Height (lowest setting) (armrest to floor)	Patient Height (highest)	Patient Height (lowest)	Weight Limit
10400-1/8	53"	45"	5'10"	5'2"	350 lbs.
10401-1/8	46"	37"	5'2"	4'6"	350 lbs.
10402-1/8	61"	53"	6'6"	5'10"	350 lbs.
10416-1	40"	31"	4'6"	4'	175 lbs.
10425-1/8	50.6"	42.6"	5'10"	4'8"	250 lbs.
10426-1/8	42.13"	36.13"	5'2"	4'6"	250 lbs.
10427-1/8	60.25"	48.25"	6'9"	5'10"	250 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.