



BARIATRIC

**850 LB**

WEIGHT CAPACITY

11138-1

## Large Bariatric Drop Arm Commode

11138-1 Assembled, 1/cs  
 11138B Replacement Back Rest, 1/ea  
 11138SEAT Replacement Seat, 1/ea  
 10374BK-8 Replacement Tip, 1" Black,  
 (Has replaced 11100B), 8bx/cs

### FEATURES:

- 850 lb. Weight Capacity
- Extra Large opening in toilet seat, 8 1/2" x 12".
- Durable heavy duty, powder coated steel construction tubing is sturdy and easy to maintain.
- Easy-to-release drop arm mechanism allows for safe lateral transfer to and from commode.
- Extra large, padded, deep seat provides comfort.
- Padded armrests and backrest provide comfort and support.
- Comes complete with extra large commode bucket and cover that removes from behind.
- **Limited Lifetime Warranty.**



BARIATRIC

**1000 LB**

WEIGHT CAPACITY

11132-1

## Bariatric Extra Large Drop Arm Commode

11132-1 Assembled, 1/cs  
 11132SEAT Replacement Seat, 1/ea  
 10374BK-8 Replacement Tip, 1" Black,  
 (Has replaced 11100B), 8bx/cs

### FEATURES:

- 1000 lb. Weight Capacity
- Extra large opening in toilet seat, 8 1/2" x 12".
- Heavy duty, white powder coated steel construction.
- Extra large, padded deep seat provides comfort.
- Easy-to-release drop arm mechanism allows for safe lateral transfer to and from commode.
- Padded armrests provide comfort and support.
- Comes complete with extra large commode bucket and cover that removes from behind.
- **Limited Lifetime Warranty.**

### SPECIFICATIONS

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11138-1	25.5"	27.5"	24"	21"	18.5"	18.5"	37 lbs.	41 lbs.	850 lbs.
11132-1	36.5"	40"	36"	21"	18.5"	NA	49 lbs.	53 lbs.	1000 lbs.

commodes

drive

To order call toll free: 877.224.0946

WWW.DRIVEMEDICAL.COM