

## Folding Steel Commode


- 11148-4 Has Replaced 11145-4, 4/cs  
 11148-1 Has Replaced 11105KD-1  
 and 11145-1, 1/cs  
 11148KDR-1 Retail Box, 1/cs  
 11100G Replacement Tip, 1" Grey, 4/bx

### FEATURES:


- Easily opens and folds.
- Folds flat for convenient storage and transportation.
- Durable plastic snap on seat and lid.
- Grey powder coated steel welded construction increases strength and durability.
- Plastic arms for added comfort.
- Comes complete with 7.5 qt. commode bucket, cover and splash shield.
- **Limited Lifetime Warranty.**



**TOOL FREE  
ASSEMBLY**



**OPENS & FOLDS  
IN SECONDS  
WITH  
SNAP BUTTONS.**



**BREAKS DOWN ALMOST COMPLETELY  
FLAT FOR RETAIL BOX**



**1. PART A OPENS UP AND ASSEMBLES (SEE CIRCLES ABOVE).**



**2. PART B SNAPS INTO PART A (SEE CIRCLES ABOVE).**



**3. PART C DROPS INTO PART A**



**4. TWO SEPARATE PIECES OF PART D SNAP ONTO PART A (SEE CIRCLES ABOVE).**



### SPECIFICATIONS

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11148-4/1	18"	22.5"	13.5"	15"	22.5"	16.5"	14.8 lbs.	63 lbs.	350 lbs.
11148KDR-1	18"	22.5"	13.5"	15"	22.5"	16.5"	14.8 lbs.	16.2 lbs.	350 lbs.

**drive** commodes  
 MEDICAL DESIGN & MANUFACTURING  
 To order call toll free: 877.224.0946