



12001-4/12001KD-1

## Toilet Safety Frame

- 12001-4 Assembled, 4/cs  
 12001KD-1 Knock Down, 1/cs  
 12000-4 Knock Down, Plastic Arms, 1 per Retail Box/4 Boxes Per Case  
 12001BR-1 Toilet Safety Frame Bracket, 1/ea

### FEATURES:

- Anodized aluminum is sturdy and lightweight.
- Waterfall armrests provide additional comfort and support.
- Powder coated aluminum bracket easily attaches frame to toilet.
- Arms are height and width adjustable.
- Limited Lifetime Warranty.



**THIS ITEM IS AVAILABLE IN A 1-PACK!**

**COMES WITH ADJUSTABLE (2 POSITION) WIDTH PADDED ARMS!**



**FOR USE WITH 12001-4 AND 12000-4 TOILET SAFETY FRAMES.**



12000-4

## Plastic Tub Rail

- 12039 12" White, 1/bx  
 12044 16 1/2" White, 1/bx

### FEATURES:

- Attractive, ergonomically designed tub rail provides a comfortable, safe grip.
- Tool free, easy-to-install.
- Width adjusts from 2 3/4" - 6 1/2".
- Attractive retail packaging.
- Limited Lifetime Warranty.



**TOOL FREE INSTALLATION**



12044



12039

### SPECIFICATIONS

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

Item #	Width Between Arms (smallest)	Width Between Arms (largest)	Height Setting (highest)	Height Setting (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12001-4	18.5"	22.75"	32"	28"	4 lbs.	18 lbs.	300 lbs.
12000-4	19"	23.50"	32"	28"	4 lbs.	20 lbs.	300 lbs.
12001KD-1	18"	25"	32"	28"	4 lbs.	5 lbs.	300 lbs.

Item #	Adjustable Clamp (smallest)	Adjustable Clamp (widest)	Product Height	Product Weight	Carton Shipping Weight	Weight Limit
12039	2.75"	6.5"	12"	3.6 lbs.	5.5 lbs.	300 lbs.
12044	2.75"	6.5"	16.5"	3.9 lbs.	5.8 lbs.	300 lbs.

bath safety

To order call toll free: 877.224.0946

**drive**  
 WWW.DRIVEMEDICAL.COM