

## New & Improved

### K.D. Padded Transfer Bench

- 12005KD-2 Knock Down Padded, Tool Free Assembly Back, Legs and Arm, (Has replaced 12005-2), 2/cs  
 12005KD-1 Knock Down Padded, Tool Free Assembly Back, Legs and Arm, (Has replaced 12005-1), 1/cs  
 12005KDR-1 Knock Down Padded, Tool Free Assembly Back, Legs and Arm, Retail, 1/cs  
 12011SS Replacement Suction Cup for Transfer Bench, Small, 1pr/bx  
 12011SL Replacement Suction Cup with Screws for Transfer Bench, Large, 1pr/bx  
 12005TP Replacement Transfer Pad, 1/ea

#### FEATURES:

- Comfortable cushioned seat and backrest.
- New "A" frame construction provides additional stability.**
- 1" Aluminum frame is lightweight, sturdy and corrosion resistant.
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs.
- Back reverses without tools.
- Extra large suction cups provide added safety.
- Limited Lifetime Warranty.



**MUST PULL "TAB" ON SUCTION CUP TO RELEASE FROM THE TUB.**

12005KD-2

**SUPPORT COLLAR REDUCES RATTLE!**



**ASSEMBLED WITH TOOL FREE REMOVABLE BACK**



#### SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12005KD-2	27.5"	24"	16.5"	16"	21.75"	17.75"	11.68 lbs.	26.3 lbs.	400 lbs.
12005KD-1*	27.5"	24"	16.5"	16"	21.75"	17.75"	11.68 lbs.	13.9 lbs.	400 lbs.
12005KDR-1	27.5"	24"	16.5"	16"	21.75"	17.75"	11.68 lbs.	15.8 lbs.	400 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

bath safety

To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM