




LEGS AND BACK EASILY ASSEMBLED WITHOUT TOOLS!



SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12011KD-2	33"	26"	24"	18.5" - 19.5"	21.75"	17.75"	10.3 lbs.	23.8 lbs.	400 lbs.
12011KD-1	33"	26"	24"	18.5" - 19.5"	21.75"	17.75"	10.30 lbs.	12.8 lbs.	400 lbs.
12011KDR-1	33"	26"	24"	18.5" - 19.5"	21.75"	17.75"	10.30 lbs.	13.8 lbs.	400 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.



ENVIRONMENT FRIENDLY PRODUCT
 Recycled materials were used to make this product.

New & Improved

K.D. Plastic Transfer Bench

- 12011KD-2** Knock Down Tool Free Back, Legs and Arm, 2/cs (Has replaced 12011-2 and 12020-2)
- 12011KD-1** Knock Down Tool Free Back, Legs and Arm, 1/cs (Has replaced 12011-1)
- 12011KDR-1** Knocked Down, Retail, Tool Free Back, Legs and Arm, 1/cs
- 12011SS** Replacement Suction Cup for Transfer Bench, Small, 1pr/bx
- 12011SL** Replacement Suction Cup with Screws for Transfer Bench, Large, 1pr/bx
- 12011-A** Replacement Arm for 12011KD and 12011 1/ea

FEATURES:

- New "A" frame construction provides additional stability.
- Durable blow-molded plastic bench and backrest.
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs.
- All transfer benches are reversible to accommodate any bathroom.
- No exposed hardware to injure patient.
- Extra large, locking suction cups provide added safety.
- Tool-Free Assembly requires connecting Back, Legs and Arm.
- Limited Lifetime Warranty.