



ASSEMBLED WITH TOOL
FREE REMOVABLE BACK



ENVIRONMENT

FRIENDLY PRODUCT

Recycled materials were
used to make this product.



12490-2*



12029-1



ASSEMBLED
WITH
TOOL FREE
REMOVABLE
BACK

12021C



EURO-STYLE
CLIP FOR
HEIGHT
ADJUSTMENT!



12021-4

"U" Shaped Shower Chair

12490-2*

With Back, Assembled with Tool Free Removable Back,
(Has replaced 12490KD-6), 2/cs
Replacement Tips for 12490, 4/bx

12002T

FEATURES:

- Contoured, "U" Shape design makes bathing more comfortable.
- Easy to install, Tool Free assembly of back support.
- Back can be adjusted to 1 of 3 positions for maximum comfort.
- Blow molded construction provides strength and durability.
- Aluminum frame is lightweight and corrosion proof.
- Angled legs provide additional stability and safety.
- Large suction style cups reduce slipping.
- Limited Lifetime Warranty.

Deluxe Bariatric Bath Bench with Cross Frame Brace

12028-2*

Without Back, Fully Assembled, 2/cs

12029-1

With Tool Free Removable Back, Fully Assembled, 1/cs

12029-2

With Tool Free Removable Back, Fully Assembled, 2/cs

12028T

Replacement Tips for 12028, 12029, 4/bx

FEATURES:

- 600 lb. Weight Capacity.
- Tool free removable back.
- Manufactured with 1" aluminum tubing.
- Frame has width/depth cross braces attached with aircraft-type rivets.
- Large 16" seat depth.
- Durable, blow molded plastic bench.
- Drainage holes reduce slipping.
- Limited Lifetime Warranty.



THIS ITEM
IS AVAILABLE
IN A 1-PACK!

Deluxe Bariatric Bath Bench

12021-4

With Back, Assembled, 4/cs

12021KD-1

With Back, Knock Down, 1/cs

12022-4

Without Back, Assembled, 4/cs

12022KD-1

Without Back, Knock Down, 1/cs

12021C

Replacement Euro-Style Clip, 1/ea

12021T

Replacement Tips for 12021, 4/bx



THIS ITEM
IS AVAILABLE
IN A 1-PACK!

FEATURES:

- 500 lb. Weight Capacity.
- Cross brace attached with aircraft type rivets.
- Blow molded bench and back provides comfort and strength.
- Drainage holes in seat reduce slipping.
- Adjustable height legs.
- Aluminum frame is lightweight, durable and corrosion proof.
- Limited Lifetime Warranty.

SPECIFICATIONS

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12490-2	18"	16"	22.5"	15"	17.75"	13.75"	6.0 lbs.	13.5 lbs.	300 lbs.
12028-2	22"	17.25"	18.75"	16"	20"	16"	8.5 lbs.	19 lbs.	600 lbs.
12029-1/2	22"	17.25"	18.75"	16"	20"	16"	10 lbs.	11/21.5 lbs.	600 lbs.
12021-4	18"	20"	17.5"	12"	21"	14"	7.4 lbs.	32.6 lbs.	500 lbs.
12021KD-1	18"	20"	17.5"	12"	21"	14"	7.4 lbs.	8.6 lbs.	500 lbs.
12022-4	18"	20"	17.5"	12"	21"	14"	6.0 lbs.	27 lbs.	500 lbs.
12022KD-1	18"	20"	17.5"	12"	21"	14"	6.0 lbs.	8.0 lbs.	500 lbs.