



**TOOL FREE  
REMOVABLE  
BACK, LEGS  
& ARM!**



**THIS ITEM  
IS AVAILABLE  
IN A 1-PACK!**



**ENVIRONMENT FRIENDLY PRODUCT**

Recycled materials were used  
to make this product.

12021C



**EURO-STYLE  
CLIP ENSURES  
A SECURE FIT!**



12025KD-2



**MUST PULL "TAB"  
ON SUCTION CUP  
TO RELEASE FROM  
THE TUB.**

**BARIATRIC  
500 LB  
WEIGHT CAPACITY**

## K. D. Bariatric Transfer Bench

- 12025KD-2 Tool Free Assembly,  
(Has replaced 12025-2), 2/cs
- 12025KD-1 (Has replaced 12025-1  
and 12009-1), 1/cs
- 12011SS Replacement Suction Cup for  
Transfer Bench, Small, 1pr/bx
- 12011SL Replacement Suction Cup for  
Transfer Bench, Large, 1pr/bx

### FEATURES:

- 500 lb. Weight Capacity.
- Tool Free Assembly.
- Durable blow-molded plastic bench and backrest.
- 1" Aluminum frame is lightweight, sturdy and corrosion resistant.
- Height adjusts in 1/2" increments with unique "Dual Column" extension legs.
- Backrest is reversible to accommodate any bathroom.
- Extra large locking suction cups provide added safety.
- Limited Lifetime Warranty.



12009-2

**EXTRA MIDDLE LEGS  
FOR ADDED STRENGTH.**

**BARIATRIC  
600 LB  
WEIGHT CAPACITY**

## Deluxe Bariatric Transfer Bench with Cross Frame Brace

- 12009-2 Fully Assembled, 2/cs
- 12011SS Replacement Suction Cups  
for 12009, 6/bx

### FEATURES:

- 600 lb. Weight Capacity.
- Manufactured with 1" aluminum tubing.
- Durable, blow molded plastic bench and backrest.
- Limited Lifetime Warranty.

### SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12025KD-1/2	33"	26.5"	25"	16.5"	23.25"	18.25"	12 lbs.	13.8/27 lbs.	500 lbs.
12009-2	32"	27"	22.25"	16"	21"	16"	14 lbs.	32 lbs.	600 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

bath safety

To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM