



12030



LOCKING
KNOB



12038-3

Open, Padded Raised Toilet Seat with Four Locking Brackets

12030
12030B

1/cs (Has replaced 12031 with plastic seat)
Replacement Bracket for 12030, 4/bx

FEATURES:

- Padded seat provides extra comfort.
- Open seat promotes patient independence.
- Four adjustable locking clamps provides versatility.
- Adjustable angle and seat height.
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.

**RAISED SEAT ADJUSTS TO DIFFERENT ANGLES
AND THE OPENING CAN BE POSITIONED AT
90° INCREMENTS TO THE FRONT**

New & Improved Raised Toilet Seat

12038-3
12040-3

Bulk Pack, 1/polybag, 3/cs
Retail, 3 retail boxes/cs, 3/cs

FEATURES:

- Designed for individuals who have difficulty sitting down, or getting up from the toilet.
- Composite material is sturdy, safe and easy-to-clean.
- Color: White.
- Fits most toilets.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Width (seat)	Depth (seat)	Seat Height	Product Weight	Carton Shipping Weight	Weight Limit
12030	16"	15.5"	5"	6 lbs.	9 lbs.	300 lbs.
12038-3	15"	15"	6"	2.4 lbs.	8.7 lbs.	300 lbs.
12040-3	15"	15"	6"	2.4 lbs.	10.7 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.