



**THIS ITEM  
IS AVAILABLE  
IN A 1-PACK!**



## ENVIRONMENT FRIENDLY PRODUCT

Recycled materials were used  
to make this product.



### Deluxe, K.D. Aluminum Bath Seat

- 12202KD-1\* With Tool Free Removable Back, 1/cs  
12202KD-4 With Tool Free Removable Back  
(Has replaced 12202-4, 12105-4 and 12002KD-6), 4/cs  
12202KDR-1 With Tool Free Removable Back, Retail,  
(Has replaced 12105KDR-1, 12002KDR-1, 12002KDR-2,  
12202KD-2 and 12495KD-1), 1/cs  
12202KDRT-1 With Tool Free Removable Back,  
Retail, Teal, (Has replaced 12002T-2), 1/cs  
12202KDRY-1 With Tool Free Removable Back,  
Retail, Yellow, (Has replaced 12002Y-2), 1/cs  
12202KDRB-1 With Tool Free Removable Back,  
Retail, Blue, (Has replaced 12002B-2), 1/cs  
12002T Replacement Tips for 12202, 4/bx

#### FEATURES:

- Easy Snap together assembly, back, seat and legs.
- Drainage holes in seat reduces slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Angled legs with suction style tips provide additional stability.
- Limited Lifetime Warranty.

12202KD-4



**SUPPORT  
COLLAR  
PREVENTS LEG  
MOVEMENT!**



12202KDRT-1



12202KDRY-1



12202KDRB-1



**LEGS AND  
BACK EASILY  
ASSEMBLED  
WITHOUT  
TOOLS!**

### Bath Bench Arms

12002A\* 1/pr

#### FEATURES:

- Easily attaches to removable back bath bench to provide additional stability. No tools required.
- For use with 12202-4 bath bench.



### Deluxe Handheld Shower Spray and Diverter Valve

12037 White, with Diverter Valve, 1/bx  
12046 Diverter Valve, (sold separately) 1/ea

#### FEATURES:

- On/Off switch built into handle for easy access.
- Wall holder and extra-long 80" white reinforced nylon hose provides extra convenience.
- Comes with diverter valve that allows either the regular shower head to operate or be used as a hand held shower spray.
- Attractive retail packaging.
- Limited Lifetime Warranty.



#### SPECIFICATIONS

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12202KD-1/4	17"	20"	17"	14.5"	19.75"	15.75"	6.5 lbs.	7.5/29 lbs.	400 lbs.
12202KDR B, T, Y-1	17"	20"	17"	14.5"	19.75"	15.75"	6.5 lbs.	8.5 lbs.	400 lbs.

bath safety

To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM