



HCPCS
E0143/E0156



Two Wheeled Walker with Seat

2 Button

| Item # | Color | HCPCS | UOM |
|--------|-------|---------------|------|
| 1239BL | Blue | E0143 / E0156 | 1/ea |
| 1239RD | Red | E0143 / E0156 | 1/ea |

- Features comfortable and convenient padded seat and backrest.
- Quickly folds when not in use.
- Rear leg tips act as brakes when pressed down.
- Height adjustable handles.
- Large swivel wheels for added maneuverability.
- Contoured handles for a comfortable grip.

specifications

| | |
|------------------------|----------------------|
| Overall Width | 20.75" |
| Depth | 23.25" |
| Handle Height | 29" - 38" |
| Seat Height | 20.5" |
| Seat Dimensions | 15.875" (W) x 8" (D) |
| Weight | 12 lbs. |
| Carton Shipping Weight | 14 lbs. |
| Weight Limit | 300 lbs. |
| Warranty | Limited Lifetime |