

**18" WIDTH
BETWEEN ARMS!**



12013

**FITS ON
ELONGATED
TOILETS!**



**NOW WITH TOOL FREE
REMOVABLE ARMS!**



**TOOL FREE
REMOVABLE
& ARMS!**

**19" WIDTH
BETWEEN ARMS**



New & Improved

Premium Seat Rizer with Removable Arms

12402

Standard Toilets (Has replaced 12400
and 12403), 1/cs

FEATURES:

- Allows individual to use their existing toilet seat and lid.
- Tool free removable arms offers flexibility to the users.
- Blow molded construction offers a durable contemporary look.
- Easily attaches to toilet.
- Color: White.
- Limited Lifetime Warranty.



SPECIFICATIONS

Item #	Width (seat)	Depth (seat)	Seat Height	Width Between Arms	Product Weight	Carton Shipping Weight	Weight Limit
12013	17"	16.5"	4.5"	18"	4.4 lbs.	6.5 lbs.	300 lbs.
12402	14"	19"	3.5"	19"	3.4 lbs.	4.4 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.