



All Purpose Stool with Adjustable Arms

12455

White, 1/cs



FEATURES:

- Stool features a padded seat and back for added comfort.
- Angled seat makes sitting down and getting up easy.
- Manufactured with white powder coated steel for an attractive, durable finish. Legs are chrome plated steel.
- Comes with removable arm supports with adjustable width for added safety and comfort.
- Easily assembled.
- Adjustable height to accommodate most individuals.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Width (outside legs)	Width (seat)	Depth (outside legs)	Depth (seat)	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12455	17"	14"	17"	10.5"	27"	21"	12.5 lbs.	15.5 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.