



12004KD-1



**THIS ITEM
IS AVAILABLE
IN A 1-PACK!**



**SUPPORT COLLAR
REDUCES RATTLE!**



12004KDRY-1



12004KDRT-1



12004KDRB-1



K.D. Bath Stool

- 12004-2* White, Assembled, 2/cs
12004KD-1 Knocked Down, White,
(Has replaced 12057, 12058, 12059), 1/cs
12004KDRB-1 Knocked Down Retail, Blue,
(Has replaced 12004B-2), 1/cs
12004KDRT-1 Knocked Down Retail, Teal,
(Has replaced 12004T-2), 1/cs
12004Y-2* Yellow, 2/cs
12004KDRY-1 Knocked Down Retail, Yellow, 1/cs
12002T Replacement Tips for 12204, 4/bx

FEATURES:

- Manufactured with 1" aluminum tubing to ensure lightweight stability.
- Impact resistant, composite seat is crack proof and tarnish resistant.
- Limited Lifetime Warranty.



12487



12486

Folding Bath Bench

- 12486 Without Back, Retail, 1/cs
12487 With Back, Retail, 1/cs
12486T Replacement Tips for 12486, 12487, 4/bx

FEATURES:

- Blow molded bench provides comfort and strength.
- Drainage holes in seat reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Angled legs provide additional stability.
- Limited Lifetime Warranty.



**CONVENIENTLY
FOLDS FLAT!**

SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
12004-2/Y	13.5"	13"	13.5"	13"	17"	13"	3.5 lbs.	8.2 lbs.	300 lbs.
12004KDRB, T, Y-1	13.5"	13"	13.5"	13"	17"	13.5"	3.5 lbs.	5.5 lbs.	300 lbs.
12486	17.5"	20"	20"	12"	18"	16"	4.4 lbs.	6.6 lbs.	300 lbs.
12487	17.5"	20"	20"	12"	18"	16"	5.0 lbs.	7.2 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

bath safety

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM