



## Pivot & Tilt Overbed Table

13000

1/cs

### FEATURES:

- The mast of the table pivots and can be locked in one of the three positions from flat on the floor to 90°.
- The pivot feature allows the table top to be positioned closer to individuals in a bed, sitting in a wheelchair, or folded for storage.
- The table top tilts with the angle of the mast to provide a flat surface regardless of the angle of the mast.
- Four 2" casters (2 locking).
- Color: grey table top, black frame.
- Limited Lifetime Warranty.

### NEW & IMPROVED PACKAGING!



**ONLY WHEELS HAVE  
TO BE INSTALLED.**

### SPECIFICATIONS

Item #	Tabletop Adjustment (highest)	Tabletop Adjustment (lowest)	Tabletop Surface (width)	Tabletop Surface (depth)	Overall Width at Base	Overall Depth at Base	Product Weight	Carton Shipping Weight	Weight Limit
13000	37"	28"	23.5"	15.75"	25"	25"	21 lbs.	24 lbs.	40 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.