



SILVER VEIN FINISH ON 13067

EASY TO MAINTAIN &
LESS INSTITUTIONAL
LOOKING!



TOP TILTS 33°
IN EITHER
DIRECTION

IDEAL FOR HOMECARE!

Non-Tilt Overbed Table



- 13003 Chrome Base and Mast, Walnut Top, 1/cs
13067 Silver Vein Base and Mast, Walnut Top,
(Has replaced 13007), 1/cs

FEATURES:

- Walnut, wood grain low pressure laminate top.
- Table top can be raised or lowered in infinite positions between 28"-45".
- Chrome-plated steel "H" base provides security and stability. (Model # 13003)
- Silver Vein "H" Base provides security and stability. (Model # 13067)
- 15" x 30" table top is locked securely when height adjustment handle is released.
- Table top can be raised with slightest upward pressure.
- 2" swivel casters allow for easy maneuverability.
- Limited Lifetime Warranty.

ASK ABOUT "U" BASE OPTION

Deluxe, Tilt-Top Overbed Table

- 13008 1/cs
13003NC Caster for Overbed Table, without Lock,
1/ea (For use with 13003, 13007, 13008,
13067BV, 13064 Overbed Table)
13003NCL Caster for Overbed Table, with Lock,
1/ea (For use with 13064 overbed tables)

FEATURES:

- Easily assembled.
- Walnut, wood grain low pressure laminate top.
- Top can be tilted 33° in either direction.
- Top can be raised or lowered in infinite settings between 26.5"-39".
- Swivel casters allow for easy maneuverability.
- Chrome-plated steel "H" base provides security and stability.
- Table top is locked securely when height adjustment handle is released.
- Table top can be raised with the slightest upward pressure.
- Limited Lifetime Warranty.

SPECIFICATIONS

Item #	Tabletop Adjustment (highest)	Tabletop Adjustment (lowest)	Tabletop Surface (width)	Tabletop Surface (depth)	Overall Width at Base	Overall Depth at Base	Product Weight	Carton Shipping Weight	Weight Limit
13003	45"	28"	30"	15"	26.5"	15.5"	22 lbs.	25 lbs.	40 lbs.
13008	46.25"	30.5"	30"	15"	26.5"	15.5"	25 lbs.	28 lbs.	50 lbs.
13067	45"	28"	30"	15"	26.5"	15.5"	22 lbs.	25 lbs.	40 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.