



**SILVER VEIN FINISH AVAILABLE**  
EASY TO MAINTAIN &  
LESS INSTITUTIONAL  
LOOKING!



**TOP TILTS 33°  
IN EITHER  
DIRECTION**

## IDEAL FOR HOMECARE!

### Non-Tilt Overbed Table

**13003** Chrome Base and Mast, Walnut Top, 1/cs  
**13067** Silver Vein Base and Mast, Walnut Top,  
(Has replaced 13007), 1/cs

#### FEATURES:

- Walnut, wood grain low pressure laminate top.
- Table top can be raised or lowered in infinite positions between 28"-45".
- Chrome-plated steel "H" base provides security and stability.
- 15" x 30" table top is locked securely when height adjustment handle is released.
- Table top can be raised with slightest upward pressure.
- Base and mast are available in chrome or silver vein finishes.
- 2" swivel casters allow for easy maneuverability.
- Limited Lifetime Warranty.

**ASK ABOUT "U" BASE OPTION**

### Deluxe, Tilt-Top Overbed Table

**13008** 1/cs  
**13003NC** Caster for Overbed Table, without Lock, 1/ea (For use with 13003, 13007, 13008, 13067BV, 13064 Overbed Table)  
**13003NCL** Caster for Overbed Table, with Lock, 1/ea (For use with 13064 overbed tables)

#### FEATURES:

- Easily assembled.
- Walnut, wood grain low pressure laminate top.
- Top can be tilted 33° in either direction.
- Top can be raised or lowered in infinite settings between 26.5"-39".
- Swivel casters allow for easy maneuverability.
- Chrome-plated steel "H" base provides security and stability.
- Table top is locked securely when height adjustment handle is released.
- Table top can be raised with the slightest upward pressure.
- Limited Lifetime Warranty.

#### SPECIFICATIONS

Item #	Tabletop Adjustment (highest)	Tabletop Adjustment (lowest)	Tabletop Surface (width)	Tabletop Surface (depth)	Overall Width at Base	Overall Depth at Base	Product Weight	Carton Shipping Weight	Weight Limit
<b>13003</b>	45"	28"	30"	15"	26.5"	15.5"	22 lbs.	25 lbs.	40 lbs.
<b>13008</b>	46.25"	30.5"	30"	15"	26.5"	15.5"	25 lbs.	28 lbs.	50 lbs.
<b>13067</b>	45"	28"	30"	15"	26.5"	15.5"	22 lbs.	25 lbs.	40 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.