


A



A

## Trapeze Bar with Brown-Vein Finish

Item #	Product Description	UOM
13009KT-BV	Includes Trapeze Bar, Triangle Chains and Clamps	1/cs

- Heavy-gauge two-piece square steel tubing provides maximum support.
- Mounts easily to metal frame headboards and fits most major manufacturers'.
- Attractive brown-vein finish. (Figure A)
- Easy to install without tools.
- Easy-to-adjust height, length and horizontal position.
- Can be shipped UPS.

### specifications

User Weight Limit	250 lbs.
Warranty	Limited Lifetime

### accessories

13009BV	Bar Only	1/cs
13009CBV	Clamp for Trapeze, 2 Pieces	1/cs

## Trapeze Base with Brown-Vein Finish

Item #	Product Description	UOM
13017BV	For Use by Patients up to 250 lbs.	1/cs

- For use when metal headboards are not available.
- Attractive brown-vein finish. (Figure A)
- Folds easily for storage and transport.
- Easy to assemble. No tools necessary.
- Adjustment levers on legs.
- For use with home-style beds.

### specifications

Weight	24 lbs.
Carton Shipping Weight	27.5 lbs.
User Weight Limit	250 lbs.
Warranty	Limited Lifetime