



Economy I.V. Pole, 4 Leg with Removable Top

13029	4 Hook, 1/cs
13033	2 Hook, 1/cs (Has replaced 13015)
13033SV	Silver Vein, 2 Hook, 1/cs

COMPARE TO:	
drive	Brewer®
13029	11300

FEATURES:

- 4 rubber wheels provide for a smooth transport and maneuverability.
- Removable hooks make this economy I.V. pole versatile.
- Easily converted to 2 hook or 4 hook with easy to release push-pin.
- Chrome plated steel with weighted base provides strength, durability and reduces risk of tipping.
- Height range of 40"-82".
- Locking collar allows for easy height adjustment.
- Limited Lifetime Warranty.




Deluxe Heavyweight I.V. Pole, 5-Leg, Infusion Pump Style

13022	1/cs
-------	------

COMPARE TO:		
drive	Brewer®	Clinton®
13022	43416	IV-33

FEATURES:

- Heavy center weight provides added stability.
- Ideal for use with multiple pumps.
- 5-leg base with ball bearing casters provide optimum versatility.
- Chrome plated pole with easy to adjust lock allows height of pole to be varied.
- Comes with 2 hooks.
- Limited Lifetime Warranty.



I.V. Pole Top (Only)

13028	4 Hook, 1/cs (shown)
13040	2 Hook, 1/cs

FEATURES:

- Easily attaches to 13016, 13022, 13029, 13033 and 13043 I.V. poles.
- Chrome plated steel is easy to maintain.
- Limited Lifetime Warranty.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

patient room

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM