

13030-1SV



**SILVER VEIN FINISH**  
EASY TO MAINTAIN &  
LESS INSTITUTIONAL LOOKING!

COMPARE TO:  
**drive** Brewer® Clinton®  
 13030-1SV 11200-1 T-42



THIS ITEM  
IS AVAILABLE  
IN A 1-PACK!

## Foot Stool

13030-1SV (Has replaced 13030-1), 1/cs  
 13030-2SV (Has replaced 13030-2), 2/cs

### FEATURES:

- Steel welded silver vein construction.
- Non-slip ribbed rubber platform.
- Reinforced rubber tips.
- Limited Lifetime Warranty.

COMPARE TO:  
**drive** TFI® Brewer®  
 13037-1SV 6354 T-42

**BARIATRIC**  
**500 LB**  
 WEIGHT CAPACITY

## Bariatric Foot Stool

13037-1SV (Has replaced 13037-1), 1/cs  
 13062-1SV With Handrail,  
 (Has replaced 13062-1), 1/cs

### FEATURES:

- 500 lb. weight capacity.
- Manufactured with 1" steel tubing.
- Non-skid ribbed rubber platform.
- Cross brace attached with aircraft-type rivets.
- Limited Lifetime Warranty.



**SILVER VEIN FINISH**  
EASY TO MAINTAIN & LESS  
INSTITUTIONAL LOOKING!

**NEW  
KNOCK  
DOWN  
DESIGN!**




13031-1SV

COMPARE TO:  
**drive** TFI® Brewer®  
 13031-1SV 6351 11220



THIS ITEM  
IS AVAILABLE  
IN A 1-PACK!

**SILVER VEIN FINISH**  
EASY TO MAINTAIN &  
LESS INSTITUTIONAL  
LOOKING!



## Foot Stool with Handrail

13031-1SV Deluxe with handrail,  
 Assembled  
 (Has replaced 13031-1), 1/cs  
 13031-2SV Deluxe with handrail,  
 Assembled  
 (Has replaced 13031-2), 2/cs

### FEATURES:

- Handrail provides additional support and safety.
- Steel welded construction.
- Non-slip ribbed rubber platform.
- Reinforced rubber tips.
- Limited Lifetime Warranty.

## SPECIFICATIONS

Item #	Step Surface (width)	Step Surface (depth)	Step Height	Overall Width	Overall Depth	Handle Height	Product Weight	Carton Shipping Weight	Weight Limit
13030-1SV	14"	11"	9"	17"	13"	NA	5 lbs.	7 lbs.	300 lbs.
13030-2SV	14"	11"	9"	17"	13"	NA	5 lbs.	13 lbs.	300 lbs.
13037-1SV	14"	10"	9"	18.5"	14"	NA	7 lbs.	9 lbs.	500 lbs.
13062-1SV	14"	10"	9"	18.5"	14"	34"	10 lbs.	13 lbs.	500 lbs.
13031-1SV	14"	11"	9"	17"	13"	34"	8 lbs.	11 lbs.	300 lbs.
13031-2SV	14"	11"	9"	17"	13"	34"	8 lbs.	19 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.