



INSTALLS ON  
VIRTUALLY ALL  
SMOOTH SURFACES  
WITHOUT TOOLS!



## Deluxe, International Grade, Adjustable Length Suction Cup Grab Bars

- 13063 Single Hand Grab Bar with Suction Cups, 1/bx
- 13063S Adj. Grab Bar with Suction Cups, Small, 1/bx
- 13063M Adj. Grab Bar with Suction Cups, Medium, 1/bx
- 13063L Adj. Grab Bar with Suction Cups, Large, 1/bx
- 13063P Replacement Suction Cup End, 1/ea

### FEATURES:


- Large rubber suction cups (4.5") provide superior adhesion properties.
- Portable, quickly and easily installed on virtually any smooth surface without tools. No holes to drill.
- Provides assistance and added safety to the home or when traveling.
- For use on a shower wall, in the bath mounted on the tub, for use as a temporary towel bar, near a chair or stairwell for assistance when standing or climbing stairs.
- Telescoping bar for one or two handed use.
- Available in one fixed length and three telescoping lengths.
- Stainless steel release and locking levers make installing and removing the suction cup grab bars easy.
- Color: White
- Warranty: 3 years on the grab bar, 1 year on the suction cups.



13063M



**ORDER**  
**ITEM # PLANOGRAMB**  
**CONTAINS**  
(4) 13063, (8) 13063S,  
(4) 13063M, (4) 13063L  
AND (1) HEADER SIGN




**CAUTION: THE SUCTION CUP  
GRAB BARS ARE TO BE  
USED AS BALANCE ASSIST  
ONLY. THEY ARE NOT FOR  
BODY WEIGHT LEVERAGE.**

## Suction Cup Grab Bar Angle Adapter

13063-ADP 2/bx

### FEATURES:

- Designed so that the grab bar can attach to inside corners or curves.
- The aluminum adapter swivels 360°.
- Easily attaches to grab bar without tools.



### SPECIFICATIONS

Item #	Length Adj. (min)	Length Adj. (max)	Product Weight	Shipping Weight
13063	12.75"	12.75"	.4 lbs	.9 lbs
13063S	17.25"	22.2"	.8 lbs	1.3 lbs
13063M	21.75"	26.75"	1.2 lbs	1.7 lbs
13063L	25.75"	30.75"	1.6 lbs	2.1 lbs

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.