



13246

REPLACEMENT PARTS:

**External Charger
with Power Cord**

13255 1/ea

**Rechargeable
Battery**

13258 1/ea

**Wall Mounted
Bracket**

13260 1/ea

Stand-Assist Lift

13246 1/cs

FEATURES:

- Provides stable assistance in standing, transferring and toileting.
- Four sling attachment points enable the lift to accommodate a variety of stand-up sling designs.
- Extra-wide, non-skid foot plate provides a stable base for foot positioning.
- **Soft, padded adjustable leg support - includes a support strap to secure legs to the support.**
- Easy roll 4" casters mean stability and easy maneuverability over a variety of surfaces.
- 4.5 under bed clearance provides access undermost long term care beds.
- 24 V DC motor provides power and reliability.
- Motor has an emergency stop button and manual lowering feature.
- Fast charging, re-chargeable battery is removable to eliminate lift down-time.
- Battery has a charged and low-charge indicator. Audible tone when charge is complete.
- Battery charge meter indicates the battery power status while in use. Audible tone when battery power is below 50%.
- Includes removable, rechargeable battery (item #13258), charging station (#13255) and wall mount bracket (item #13260).

To prevent lift downtime - order an additional re-chargeable battery, item # 13258. One for the lift and one for the charging station.

**COMPARE TO:**

Invacare®	RPS350-1
Invacare®	RPS440-1
Graham Field	LF2020

SPECIFICATIONS

Item #	Boom Maximum Height	Boom Minimum Height	Base Width Open	Base Width Closed	Base Clearance	Casters	Product Weight	Carton Shipping Weight	Weight Limit
13246	81.1"	26.2"	40.2"	24.5"	4.5"	3" (Front) 4" (Rear)	113 lbs.	128 lbs.	400 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

patient room

drive
 WWW.DRIVEMEDICAL.COM

To order call toll free: 877.224.0946