



13606

**NEW DESIGN ALLOWS FOR  
EASY INSTALLATION &  
REMOVAL OF BATTERIES.**

### Pressure Sensitive Chair and Bed Patient Alarm

- 13605 Chair Size, (Comes with #13610 Standard Alarm and Chair Size Pad), 1/ea
- 13606 Bed Size, (Comes with #13610 standard alarm and bed size pad), 1/ea
- 13608 Chair Size, (Comes with #13611 alarm with reset button and chair size pad), 1/ea
- 13609 Bed Size, (Comes with #13611 alarm with reset button and bed size pad), 1/ea



13605

EASILY ATTACHES  
TO WHEELCHAIR.



ALARM SOUNDS WHEN PRESSURE  
IS REMOVED FROM SENSOR PAD.

### Replacement Parts

- 13610 Replacement Alarm Only For Use with 13605, 13606, 1/ea
- 13611 Replacement Alarm with Reset Button For Use with 13605, 13606, 1/ea
- 13610PL Sensor Pad Only, Bed Size, (6 Month Warranty), 1/ea
- 13610PS Sensor Pad Only, Chair, (6 Month Warranty), 1/ea

### FEATURES:

- Alerts care giver with audio alarm when patient gets out of a chair or bed.
- Pressure sensitive pad connects to the audio alarm.
- The alarm comes with a safety alert, ON/OFF switch and low battery warning.
- Comes with 4 AA batteries (Included).
- Pad size: (chair) 10" x 15".
- Pad size: (bed) 11" x 30".
- Warranty: (alarm) 6 months.

patient room

To order call toll free: 877.224.0946

**drive**  
WWW.DRIVEMEDICAL.COM

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.