order call toll free: 877.224.0946

Med - Aire Alternating Pressure Mattress Replacement System with Low Air Loss

14027 Mattress and Pump System, 1/cs

14027M Mattress Only, 1/cs 14027P Pump Only, 1/cs

Replacement Bladder for 14027 14027-CELL

and 14028N, 1/ea

14027-CELL 1 Replacement Bladder with

Laser Holes for 14027 and 14028N, 1/ea

HCPCS code:

FEATURES:

8-Inch high air cells are constructed with 40% Nylon/60% PVC for durability.

Mattress can operate in an alternating or static pressure mode.

The system provides alternating pressure and is designed to be used in the prevention, treatment and management of Pressure Ulcers.

- The "WHISPER QUIET", pump produces consistent air flow.
- Special vibration isolating features insure that the patient rests easily, without feeling the pumps vibrating action.
- LED and audible alarm provide critical performance data and alert caregiver when pressure falls below 10mmHg (inner mattress pressure).
- Easy to replace individual bladders.
- Mattress has 20 air bladders, nine (9) with laser holes to achieve low air loss function, to increase patient comfort, keep the patient cool and speed recovery.
- Quick fill adapter standard.
- Spring lined air tubes prevent kinks and disruption of air flow.
- CPR valves provide fast (20-Second) deflation.

SPECIFICATIONS:

Pump:

- Power 120 VAC.
- Frequency: 60 HZ.
- Pressure Range 30 mmHg to 80 mmHg (inner mattress pressure).
- Air volume 8-LPM (Litres per Minute).
- Cycle Time: 10 Minutes.
- Functions: Static and Dynamic Features.
- Power Cord Length: 15 Ft.
- Pump Dimensions: 11" (W) x 8" (H) x 3.5" (D).
- UL: 60601-1, 1st Edition-2006.
- Maximum Current: 0.5 Amps.
- Warranty: 12 Months Limited.

14027P

PUMP ONLY

Mattress:

- Topper Material: zippered all around and quilted urethane coated nylon fabric (2-way stretch) is; water resistant, vapor permeable, machine washable and low sheer.
- Base Material: zippered PVC mesh.
- 4-Air cells at head section remain static to enhance pillow function.
- Dimensions: 36" (W) x 80" (L) x 8" (H).
- Maximum Support Capacity: 350 lbs.
- Warranty: 12-Months Limited.



MATTRESS CAN OPERATE

PART # 14427A NOW COMES WITH "QUICK FILL VALVE" FOR **USE WITH PART # 14427, QUICK INFLATE ELECTRIC PUMP** (SOLD SEPARATELY)



14427 (Electric Pump)



EZ Lock/Quick Release Tubing Connector

- Standard with 14025, 14027, 14028N, 14029, 14029 DP, 14030, 14048, 14054 and 14060.
- Locks pump and mattress together.
 - Pressure tabs allow easy and quick disconnect.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.



