



Universal Mattress Cover, Defined Perimeter

| Item # | Product Description | UOM |
|--------|---------------------|------|
| 14333 | Fits Regular Beds | 1/ea |

- Open in mid-section for safe and easy patient egress and ingress.
- Built-in medium density foam bolsters safely and comfortably define edges of bed for patient and prevent falling.
- Easy to install and secure.
- Nylon cover can be wiped clean.
- For standard beds.

| specifications | |
|-----------------|-----------------------------------|
| Fits Mattresses | 6"-8" (H) x 34"-39" (W) x 80" (L) |
| Foam Sections | 4" (H) x 5½" (W) x 29½" (L) |
| Shipping Weight | 5 lbs. |
| Cover Material | Nylon |
| Foam Density | Medium |
| Warranty | Limited Lifetime |