


COMPARE TO:

MASON® GEL-U-SEAT™,
BLUE CHIP® GEL PRO DLX™, INVACARE®
COMFORTMATE, INVACARE® - INFINITY LO
CONTOUR & FLA® POSTURA® GELFOAM™



14886

Gel "E" 3" Wheelchair Cushions with Gel Bladder

Available
1 Pack

- 14886 18" (W) x 16" (D) x 3" (H), 275 Wt Cap., 1/cs
(Formerly 14890)
14891 16" (W) x 16" (D) x 3" (H), 275 Wt Cap., 1/cs
14894 20" (W) x 16" (D) x 3" (H),
275 Wt Cap., 1/cs
HCPCS code: E2603
14902 22" (W) x 18" (D) x 3" (H), 325 Wt Cap., 1/cs
14903 24" (W) x 18" (D) x 3" (H), 500 Wt Cap., 1/cs
HCPCS code: E2604

FEATURES:

- Pressure reducing cushion designed to assist in the prevention, treatment and management of pressure ulcers.
- Viscous Gel provides optimized pressure reduction, support and comfort.
- Fire Retardant polyurethane foam shell.
- Removable and replaceable cover has water resistant "Top" and "Non-Slip" vinyl base.
- 18 Month Limited Warranty on cover and cushion.

COMPARE TO:

MASON® - TITANIUM, JAY® CARE, PRIDE® -
SYNERGY, BLUE CHIP - AMARA 400 & INVACARE®
- INFINITY GENTLE CONTOUR



14882

Molded - Skin Protection & Positioning Wheelchair Cushions with Polymer Filled Bladders & Abductor

Available
1 Pack

- 14882* 18" (W) x 16" (D) x 3" (H), 250 Wt. Cap., 1/cs
14883* 20" (W) x 16" (D) x 3" (H), 250 Wt. Cap., 1/cs
HCPCS code: E2607
14884* 22" (W) x 18" (D) x 3" (H), 350 Wt. Cap., 1/cs
HCPCS code: E2608

FEATURES:

- Unique construction of high resilient (HR) molded foam provides optimum seat stabilization, positioning, pressure reduction and correct seating posture.
- Unique Polymer bladders beneath the Ischium for pressure reduction.
- Soft, durable, sealed inner liner provides protection from incontinence.
- Removable and replaceable cover is washable, flame retardant, reduces pressure, reduces friction & shearing, is urethane coated for water resistance and vapor permeability, is constructed with a non-slip bottom.
- 1 Year Limited Warranty on cover and cushion.

SPECIFICATIONS

Item #	Width	Depth	Height	Product Weight	Carton Shipping Weight	Weight Limit
14886	18"	16"	3"	5.7 lbs.	6.7 lbs.	275 lbs.
14891	16"	16"	3"	4.4 lbs.	5.4 lbs.	275 lbs.
14894	20"	16"	3"	4.4 lbs.	5.4 lbs.	275 lbs.
14902	22"	18"	3"	12.5 lbs.	14.0 lbs.	325 lbs.
14903	24"	18"	3"	14.0 lbs.	15.5 lbs.	500 lbs.
14882	18"	16"	3"	3.9 lbs.	4.9 lbs.	250 lbs.
14883	20"	16"	3"	4.4 lbs.	5.4 lbs.	250 lbs.
14884	22"	18"	3"	8.1 lbs.	9.6 lbs.	350 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

pressure prevention

To order call toll free: 877.224.0946

drive