



Bed Transporter



Item #	For Use With	UOM
15060	15003, 15004 and 15005, 15230, 15235, 15030 and 15033 beds	1/ea

Premium Guard Gel Overlay

Item #	Product Description	Product Size	UOM
14893	Gel Overlay	76" (L) x 34" (W) x 3.5" (H)	1/ea
14901	Bariatric Gel Overlay	76" (L) x 42" (W) x 3.5" (H)	1/ea

- New and improved gel overlay has an upgraded 300 lb. weight capacity and optimizes pressure redistribution.
- New upgrade: 15 gel bladders cover entire patient contact area.
- Full length "head-to-toe" gel bladders minimizes "bottoming out" and "pitting".
- New convoluted foam topper vastly improves patient comfort and air circulation.
- 4 elastic corner straps secures overlay to mattress.
- Urethane coated nylon cover is water resistant and vapor permeable.
- Vinyl base is stain resistant and waterproof.
- Bladders contain an aqueous non-toxic gel and polystyrene beads.
- Designated head and foot sections.
- 300 lb. weight capacity (14893), 500 lb. weight capacity (14901).

Bed Extension Kit



15005EXTKIT



15030EXTKIT

Item #	For Use With	UOM
15005EXTKIT	15003 beds	1/cs
15030EXTKIT	15030, 15033, 15230 and 15235 beds	1/cs
15005EXTKIT-L	15004 and 15005 beds	1/cs

- Adds 4" to length of bed. Easy to install.
- Easy to install.