




HCPCS
E0272



Cellulose Fiber Mattress



Item #	Product Size	UOM
15007	80" (L) x 36" (W) x 6" (H)	1/ea

- 12 oz. thermal-bonded polyester fiber.
- Two-sided.
- Inverted seam
- Fluid proof
- Non-allergenic.
- Antibacterial and antifungal cover.
- Easy cleaning with mild soap and water.

Quick 'n Easy Comfort Mattress

Item #	UOM
15076	1/ea

- Premium quality Home Care mattress offering the utmost in comfort.
- Patent pending super compression technology. Packaged in a box of just 14" x 14" x 37.5".
- Super-compact shipping carton. Easy to carry. Frees up space in delivery vehicles and warehouses. (Figure A)
- Non-slip bottom and 4 corner straps keep the mattress securely in place. (Figure B) Decompresses quickly.
- Mattress core made from select grade of polyurethane.
- Carry handles make lifting and placement easier.
- Inner core encased in fire-retardant barrier, meeting regulation 16 CFR 1633.
- Non-slip bottom and 4 corner straps keep the mattress securely in place.
- Durable nylon cover is waterproof for easy use and care.

specifications

Unit Dimensions	80" (L) x 36" (W) x 6" (H)
Weight	20 lbs.
Carton Shipping Weight	22 lbs.
Weight Capacity	275 lbs.
Warranty	Limited 1 Year

Premium Foam Mattress



Item #	Product Size	UOM
15017	80" (L) x 36" (W) x 6" (H)	1/ea

- 13 oz. thermal-bonded polyester fiber.
- Two-sided.
- Interverted seam.
- Fluid proof.
- Non allergenic.
- Antibacterial and antifungal cover.
- Easy cleaning with mild soap and water.