

## COMPARE TO:

COMPARE TO INVACARE®  
SOLACE & SPAN AMERICA GEOMAT HC

# NEW!

## Therapeutic 5 Zone Support Mattress

Available  
**1 Pack**

15019 80" (L) x 36" (W) x 6" (H),  
275 Weight Cap., 1/cs

HCPCS code: E0184

### FEATURES:

- Deluxe, horizontal, cross-cut, 5-zoned foam mattress provides comfort, support and pressure redistribution.
- 210 denier nylon cover reduces friction and shear, and is water resistant and vapor permeable.
- 3 Year Limited Warranty on foam.  
12 Month Limited Warranty on cover.



# NEW!


## Gravity 7 Long Term Care Pressure Redistribution Mattress

Available  
**1 Pack**

15770 36" x 80" x 6", 1/cs  
15777 80" (L) with Elevated Perimeter  
and Cut Out, 1/cs  
15784 36" x 84" x 6", 1/cs  
15785 84" (L) with Elevated Perimeter  
and Cut Out, 1/cs  
15876 36" x 76" x 6", 1/cs  
15877 76" (L) with Elevated Perimeter  
and Cut Out, 1/cs

### FEATURES:

- Multi-layered, high density foam with die-cut top layer provides optimal pressure redistribution.
- The ultra soft heel section is uniquely designed to provide extremely low skin interface pressure.
- The mattress cover is a 2-way stretch, water resistant, vapor permeable low shear material.
- Optional firm side rails help the patient transfer in and out of bed minimizing the risk of falling.
- Meets fire code 1633.
- 7 Year Limited Warranty.



### SPECIFICATIONS

Item #	Length	Width	Height	Product Weight	Carton Shipping Weight	Weight Limit
15019	80"	36"	6"	30 lbs.	35 lbs.	275 lbs.
15770	80"	36"	6"	30 lbs.	35 lbs.	300 lbs.
15777	80"	36"	6"	30 lbs.	35 lbs.	300 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

pressure prevention

To order call toll free: 877.224.0946

**drive**