



**EASY, TOOL FREE
ASSEMBLY!**

Home Bed Assist Rail & Folding Bed Board Combo

15062

Bed Assist Board,
Retail Packaged, 1/cs

FEATURES:

- Assistance for getting in and out of bed. Can be used on either side of the bed.
- **Tri fold bed board for easy portability and storage.**
- Removable handle.
- Easy, quick, no tool assembly. Comes in 2 pieces.
- Simply slides under the mattress for extra firm support.
- Removable powder coated steel handle enables user to remove when not in use.
- Handle with the "Mid Bar", creates grip for any height use.
- Rail dimension: 11.5" (W) x 22.5" (H) x 26.5" (D)
- Weight limit: 250 lbs.



**SIMPLY SLIDES
UNDER THE
MATTRESS!**



**EASY, TOOL FREE
ASSEMBLY!
COMES IN 2 PIECES!**



**EASY, TOOL FREE
ASSEMBLY!**



Home Bed Assist Rail

15064

Bed Assist Rail, Retail Packaged, 1/cs

FEATURES:

- Assistance for getting in and out of bed. Can be used on either side of the bed.
- Removable handle.
- Easy, quick, no tool assembly. Comes in 2 pieces.
- **Safety strap wraps around mattress or box spring to ensure a secure fit.**
- Special extra long safety strap enables user to be in a single, twin, double, queen, or king size bed.
- Made of powder coated steel for long life and easy cleaning.
- Handle with the "Mid Bar", creates easy grip for any height use.
- Rail dimension: 11.5" (W) x 20.75" (H) x 21.5" (D)
- Weight limit: 250 lbs.



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

patient room

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM