



**MORE SECURE
DESIGN
ATTACHES TO
BED FRAME
AND FLOOR
WITH FOUR
ADJUSTABLE
LEGS.**



**EASY TO
ASSEMBLE!
TOOL FREE
ASSEMBLY!
COMES IN 3
PIECES!**

**LEGS TOUCH
FLOOR FOR ADDED
STRENGTH & STABILITY!**



Home Bed Side Helper

Available
1 Pack

15065R-1

Retail, (Will become
RTL15065R in 2011), 1/cs

FEATURES:

- Attaches to Hollywood style home bed frame with four legs on the floor.
- Designed to provide assistance getting in and out of home bed.
- Easy to install, no tools required.
- Rail features four easily accessible grasping points.
- Allows for normal bedding drape.
- Adjustable to fit most beds. 11" wide frame for one or two handed use.
- Durable, solid steel tubing.
- Comes with 8" mounting bracket for a secure fit.
- Rail dimensions: 11" (W) x 37" (H) x 11.5" (D)
- Weight Limit: 250 lbs.
- Limited Lifetime Warranty.

COMPARE TO:

drive TFI®
15065R-1 8750

NEW & IMPROVED!

Home Bed Assist Handle with Legs that Reach the Floor

Available
1 Pack

15063-ADJ*

Adjustable Height Bed
Assist Handle, Adj. 13.5"-16.5",
Retail Packaged, (Will
become RTL15063-ADJ
in 2011), 1/cs (Has replaced 15063)

FEATURES:

- Assistance for getting in and out of home style bed. Can be used on either side of the bed.
- Removable handle, easy, quick, no tool assembly. Comes in 3 pieces.
- Simply slides under mattress with added feature of extra "no slip foam" on base bar.
- All steel, chrome plated, solid and easy to clean.
- Portable, easy to transport. Can be adjusted to flat position for easy storage.
- Rail dimension:
19.75" (W) x 13.5" (H) x 32.75" (D)
- Weight Limit: 250 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.