



Quick N' Easy Comfort Mattress

15076 80" (L) x 36" (W) x 6" (H), 1/cs

Available
1 Pack

HCPCS Code: E0272

FEATURES:

- Premium quality homecare mattress offering the utmost in comfort.
- Patent pending super compression technology packages mattress in a box with a footprint of just 14" x 14" x 37.5".
- Decompresses quickly.
- Mattress core made from select grade of polyurethane. Offers comfort, feel and advantages of latex rubber, without latex allergen.
- Inner core encased in fire retardant barrier, meeting regulation 16 CFR 1633.
- Non-slip bottom and 4 corner straps keep the mattress securely in place.
- Two carry handles on side make lifting and placement easier.
- Durable nylon cover is waterproof for easy use and care.
- Limited 12 month warranty.

SPECIFICATIONS

Item #	Length	Width	Height	Product Weight	Carton Shipping Weight
15076	80 "	36 "	6 "	20 lbs.	22 lbs.

pressure prevention

drive[™]

To order call toll free: 877.224.0946