

**DESIGNED TO COMPLY
AS PART OF THE FDA
HOSPITAL BED SYSTEM
DIMENSIONAL AND ASSESSMENT GUIDANCE
TO REDUCE THE
CHANCE OF ENTRAPMENT**



Half Length Bed Rail, Tool Free Adjustable Width, with Brown Vein Finish

15208BV


1 pr/cs

FEATURES:

- Half rail was designed to comply as part of the FDA Hospital Bed System Dimensional and Assessment Guidance to Reduce Entrapment.
- Tool free width adjustment allows the rail to be adjusted in 2" increments from 36" to 42".
- Attractive brown vein finish.
- Constructed of 1" steel.
- Fits most beds with steel frame.
- Ideal for protection when bed is elevated.
- Not for use with "Pan-Style" bed surfaces.
- Easy to install.
- Limited Lifetime Warranty.

CHROME FINISH
ON SECOND AND FIFTH BAR
ALLOWS FOR AN EASY, SCRATCH
PROOF MOVEMENT OF
ADJUSTMENT BRACKET.

CHROME BARS TO PREVENT
CLAMP FROM SCRATCHING.




RAIL CAN BE
ADJUSTED TO INCREASE OR
DECREASE THE GAP
BETWEEN THE RAIL AND
THE MATTRESS.

BROWN VEIN FINISH
EASY TO MAINTAIN &
LESS INSTITUTIONAL!

BUTTON ALLOWS
FOR ADJUSTMENTS
FROM 36" TO 42"
TO PREVENT GAPS
BETWEEN BED AND RAIL!

CHROME BARS TO PREVENT
CLAMP FROM SCRATCHING.



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

patient room

drive
WWW.DRIVEMEDICAL.COM

To order call toll free: 877.224.0946