




BARIATRIC
600 LB
WEIGHT CAPACITY

New & Improved Full Electric Bariatric Bed

- | | |
|----------------------|--|
| 15300 | Bariatric Bed, 42" Width, 600 lb Capacity, (Has replaced 15342S), 1/cs |
| 15300BV-1HR | Bariatric Bed with 1 pair "T" Rails, (item #15021/15021BV) (Has replaced 15300-1HR) |
| 15300BV-2HR | Bariatric Bed with 2 pair "T" Rails, (item #15021/15021BV) (Has replaced 15300-2HR) |
| 15300BV-PKG | Bariatric Bed Package (contains Bariatric Bed, 42" Bariatric Foam Mattress, item #15301, and 1 pair "T" Rails, item #15021/15021BV) (Has replaced 15300-PKG) |
| 15300BV-PKG-2 | Bariatric Bed Package (contains Bariatric Bed, 42" Bariatric Foam Mattress, item #15301, and 2 pair "T" Rails, item #15021/15021BV) (Has replaced 15300-PKG-2) |



*BED PACKAGES NOW
AVAILABLE WITH
BROWN VEIN RAILS!

FEATURES:

- Quiet, smooth operation.
- Heavy duty frame ensures strength and patient safety.
- All steel construction.
- Larger sleep surface than conventional bed.
- Split pan design with removable bed ends is easy to set-up.
- No tools necessary for set-up.
- Hand control (included) provides for multiple bed positioning for patients.
- Low air loss/alternating pressure mattress available.
- Emergency manual crank included.
- Warranty: (Original purchaser) Lifetime on welds, 5 years on frame, 1 year on all other parts and components.

Bariatric Mattress

15301 **42" x 80" x 7",**
Weight Capacity 600 lbs., 1/cs

MATTRESS FEATURES:

- Top cover: Unitek 11 Nylon.
- Bottom cover: Vinyl.
- Top layer: 2.5 lb density foam.
- Middle layer: 2.5 lb density foam.
- Base: 4 lb density foam.
- Fire retardant: Compliance with 16 CFR 1633.
- Warranty: 1 year.

SPECIFICATIONS

Item #	Bed Height (min)	Bed Height (max)	Overall Length	Overall Width	Product Weight	Carton Shipping Weight	Weight Limit
15300	18.25"	26.25"	88"	42"	276 lbs.	296 lbs.	600 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.