




BARIATRIC
1000 LB
WEIGHT CAPACITY

15303

New & Improved Full Electric Bariatric Bed

- 15303** Bariatric Bed, 54" Width, 1000 lb Capacity, 1/cs
- 15303BV-1HR** Bariatric Bed with 1 pair "T" Rail, (item #15021/15021BV) (Has replaced 15303-1HR)
- 15303BV-2HR** Bariatric Bed with 2 pair "T" Rails, (item #15021/15021BV) (Has replaced 15303-2HR)
- 15303BV-PKG** Bariatric Bed Package (contains Bariatric Bed, 54" Bariatric Foam Mattress, item #15312, and 1 pair "T" Rails (item #15021/15021BV) (Has replaced 15303-PKG)
- 15303BV-PKG-2** Bariatric Bed Package (contains Bariatric Bed, 54" Bariatric Foam Mattress, item #15312, and 2 pair "T" Rails, item #15021/15021BV) (Has replaced 15303-PKG-2)



15021BV

***BED PACKAGES NOW
AVAILABLE WITH
BROWN VEIN RAILS!**

FEATURES:

- Quiet, smooth operation.
- Heavy duty frame ensures strength and patient safety.
- All steel construction.
- Larger sleep surface than conventional bed.
- Split pan design with removable bed ends is easy to set-up.
- No tools necessary for set-up.
- Hand control (included) provides for multiple bed positioning for patients.
- Low air loss/alternating pressure mattress available.
- Emergency manual crank included.
- Warranty: (Original purchaser) Lifetime on welds, 5 years on frame, 1 year on all other parts and components.

Bariatric Mattress

- 15312** 54" x 80" x 7",
Weight Capacity 1000 lbs., 1/cs

FEATURES:

- Top cover: Unitek 11 Nylon.
- Bottom cover: Vinyl.
- Top layer: 2.5 lb density foam.
- Middle layer: 2.5 lb density foam.
- Base: 4 lb density foam.
- Fire retardant: Compliance with 16 CFR 1633.
- Warranty: 1 year.

SPECIFICATIONS

Item #	Bed Height (min)	Bed Height (max)	Overall Length	Overall Width	Product Weight	Carton Shipping Weight	Weight Limit
15303	18.25"	26.25"	88"	54"	338 lbs.	363 lbs.	1000 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.