

**NEW!**



## Gravity 7 Long Term Care Pressure Redistribution Mattress

15770	36" x 80" x 6", 1/cs
15777	80" (L) with Elevated Perimeter and "Cut Out", 1/cs
15784	36" x 84" x 6", 1/cs
15785	84" (L) with Elevated Perimeter and "Cut Out", 1/cs
15876	36" x 76" x 6", 1/cs
15877	76" (L) with Elevated Perimeter and "Cut Out", 1/cs

### FEATURES:

- Multi-layered, high density foam provides optimal pressure redistribution.
- Top layer of foam is high resilient and die-cut which provides comfort and proper pressure redistribution.
- The ultra soft heel section is uniquely designed to provide extremely low skin interface pressure.
- The mattress cover is a 2-way stretch, water resistant, vapor permeable nylon.
- The cover provides a low moisture transmission rate (MVT) which keeps the patient cool and dry.
- Optional firm side rails help the patient transfer in and out of bed minimizing the risk of falling.
- Meets fire code 1633.
- **7 Year Limited Warranty.**

### SPECIFICATIONS

Item #	Length	Width	Height	Product Weight	Carton Shipping Weight	Weight Limit
15770	80"	36"	6"	30 lbs.	35 lbs.	300 lbs.
15777	80"	36"	6"	30 lbs.	35 lbs.	300 lbs.

pressure prevention

To order call toll free: 877.224.0946

**drive**