

Full Electric LTC Low Bed

- 15803SD Spring Deck for Full Electric Low Bed, 1/cs
15000HF Head and Foot Boards for LTC and Low Beds (Sold separately), 1/cs

BARIATRIC
450 LB
WEIGHT CAPACITY



SHOWN WITH
OPTIONAL
HEAD AND
FOOT BOARDS



FEATURES:

- Height adjustable sleep surface from 7.5" to 16.5".
- 3 motors provide quiet, smooth operation.
- Separate bed height and head and foot section adjustment.
- 80" x 36" spring deck surface.
- Durable powder coated epoxy finish.
- Sturdy, attractive head and foot boards.
- UL Listed.
- 450 lb. Weight Capacity.
- Warranty: (Original purchaser) Lifetime on welds and frame, two years on electronics and 1 year on all parts and components.

Accessories

Item #	15610	Wall Bumper, 1/cs
Item #	15622	Emergency Crank, 1/cs
Item #	15623	Trapeze Adapter, 1/cs
Item #	15624	Trapeze, 1/cs
Item #	15625	Assist Bar, 1/cs
Item #	15627	Traction Frame Adapter, 1/cs
Item #	15629	Panel Mounted I.V. Holder, 1/cs
Item #	15631	Side Mattress Retainer, 1 pr/cs

DECK MEASURES
LESS THAN 7.5"
IN THE LOWEST
POSITION FROM
THE FLOOR.



BED RAILS NOT
RECOMMENDED WHEN
IN THE LOW POSITION.

COMPARE TO:
Drive Graham Field®
15803SD AF4904

ACCESSORIES FEATURED ON PAGE 336 - 346
COMPLETE BED END SELECTION ON PAGE 335

SPECIFICATIONS

Item #	Bed Height (min)	Bed Height (max)	Overall Length	Overall Width	Product Weight	Carton Shipping Weight	Weight Limit
15803SD	7.5"	16.5"	87"	36"	150 lbs.	165 lbs.	450 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

institutional beds

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM