

16500CBV

New & Improved

## Tool Free Adjustable Length Home Style Bed Rail with King Size Crossbar Extension

16500BV (Has replaced 16000 and 16500SV),

1pr/cs

16500CBV (Has replaced 16000C), Brown Vein

Cross Bars with King Size Bed Extensions for 16500BV, 1pr/cs

## **FEATURES:**

- Crossbars are easy to install and remove.
  Adjusts between 36" and 72".
- Crossbars come with extensions that will accommodate a KING SIZE BED.
- Rail length can be adjusted without tools to 37", 43" or 57" long with simple push buttons.
- Fits all home style beds.
- Easily adjusts up or down with spring-loaded release.
- Attractive brown vein finish.
- Constructed of 1" steel.
- Fits all home style beds.
- Limited Lifetime Warranty,



