



HCPCS
E0570
E0574
A7006
A7013
A7015



Hercules Beetle Portable Ultrasonic Nebulizer

Item #	Product Description	UOM
18016	Hercules Beetle Portable Ultrasonic Nebulizer	1/bx

- Ultrasonic treatment provides a deep, penetrating therapy.
- Portable nebulizer with integrated, replaceable battery pack.
- Virtually silent operation.
- Easy-to-use, ultrasonic technology uses water as the conducting medium.
- Powerful nebulization offers 0.2-0.7 ml/min rate, depending on medication viscosity.
- Retail packaged.

specifications	
Method of Nebulization	Ultrasonic Vibration
Particle Size	0.5 µm to 5 µm
Unit Dimensions	3.4" (W) x 5.6" (H) x 2.0" (D)
Weight	0.51 lbs. or 0.77 lbs. w/battery
Ultrasonic Frequency	2.5MHZ Oscillation Rate
Liter Flow	Up to 0.7 ml/min
Electrical Requirements	120 VAC 60Hz
Warranty	Limited 3 Year

accessories		UOM
18016PEDMASK	Pediatric Mask	1/ea
18016ADMASK	Adult Mask	1/ea
18016MC	Medication Cup	10pk/bag
18016F	Filter	10pk/bag
18016FC	Filter Cover	1/ea
18016M	Mouthpiece	10/bag
18016N	Nasal Piece	10/bag
18016-10	Extension Tube	5/bag
18016TOP	Top Cover	1/ea
18016BASE	Base	1/ea
18016BATT	Rechargeable Battery	1/ea
18016C	A/C Charger	1/ea
18016BAG	Carry Bag	1/ea