



Voyager™ Pro

Portable Nebulizer

Item #	Product Description	UOM
18018	Voyager™ Pro Vibrating Mesh Nebulizer	1/ea

- The Voyager™ Pro is a compact and lightweight portable nebulizer which utilizes silent vibrating mesh technology to deliver treatment times in less than 10 minutes for most medications.
- Form and function needs are met by utilizing a neb cup able to accommodate up to 8ml of medication, the ability to aerosolize a wide range of medications, an auto-off function, all packaged in an attractive ergonomic design.
- Includes a mouthpiece, pediatric mask, adult mask, four AAA batteries, an AC adapter, USB port plug, and a carry bag.

specifications

Method of Nebulization	Vibrating Mesh
Particle Size	1.0 – 5.0 µm
Unit Dimensions	2.75" (W) x 5.71" (H) x 2.28" (D)
Medication Cup Capacity	2ml – 8ml
Operating Temperature	37° – 122° F
Storage Temperature	14° – 176° F
Power Requirements	Battery: 1.5V with 4 AAA Batteries Adapter: 100V – 240V / 50~60 Hz / 0.15A
Warranty	2 Year Limited Liability