



**REPLACEMENT
FILTER**

Available
1 Pack

Checker® Nebulizer

18040-Y

18040-P

NEB KIT 500

Yellow, 1/bx

Pink, 1/bx

Nebulizer Kit, 50/cs (Contains
1 each medication, cup mouth
piece and tubing)

Replacement Filter for
18040-Y and 18040-P, 1/ea

18040F

HCPSC code: E0570

FEATURES:

- Checker® cab design appeals to children.
- Available in either yellow or pink.
- Easy to replace filter is located under the front hood.
- Electrical cord located under rear hood.
- Very quiet operation.
- Powerful, dependable piston powered pump.
- Comes complete with ampule, mouth piece, pediatric aerosol mask, tubing, two replacement filters and instructions.
- **5 Year Limited Warranty.**



**LATCH HOOK DOUBLES
AS CUP HOLDER!**

**COMES WITH PEDIATRIC
AEROSOL MASK!**



18040-Y



18040-P



SPECIFICATIONS

18040

METHOD OF NEBULIZATION:	Piston Compressor
PARTICLE SIZE:	0.5µm to 5µm
UNIT DIMENSIONS:	11.8" (W) x 7.7" (H) x 4.3" (D)
WEIGHT:	3.7 lbs.
MAXIMUM PRESSURE:	30 PSI
OPERATING PRESSURE:	16 PSI
LITER FLOW:	8 LPM
ELECTRICAL REQUIREMENTS:	120 VAC 60Hz
WARRANTY:	5 Years

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

RESPIRATORY

CHAD
drive

To order call toll free: 877.224.0946