



**ASK ABOUT OUR
PRIVATE LABEL
PROGRAM!**

**OUR LIGHTEST WEIGHT
NEBULIZER AT ONLY 3.2 LBS!**



18080



18090F

**REPLACEMENT
FILTER**

Power Neb Ultra Nebulizer

- 1P** 18080 1/cs
NEBKIT 500 Disposable Nebulizer Kit, 50/cs
(Contains 1 each medication cup,
mouth piece and tubing),
(Has replaced 18060 and 18062)
- 18090F Replacement Filter, 10 per bag

HCPCS code: E0570

FEATURES:

- One of the smallest economy neb compressors available today at only 3.2 lbs!
- Available for private label (ask your local Drive rep for details).
- Built in neb kit holder for hands free use.
- Delivers consistent particle size of <5 microns and fast treatment times.
- Comes complete with disposable neb kit.
- Limited 5 Year Warranty.



18016



**REPLACEMENT
FILTER**



**Beetle Neb
18016F**

Hercules Beetle Portable Ultrasonic Nebulizer

- 1P** 18016 1/bx
(Has replaced 18015)
- 1P** 18016 ADMASK Adult, 1/ea
- 1P** 18016 PEDMASK Pediatric, 1/ea
- 18016F Replacement Filter, 10 per bag

FEATURES:

- Ultrasonic treatment provides a deep, penetrating therapy
- **PORTABLE!** Integrated, detachable Battery pack and virtually silent operation
- Easy to Use ultrasonic technology uses water as the conducting medium
Powerful nebulization offers 0.2-0.7 ml/min rate, depending on medication viscosity
- Limited 3 Year Warranty

SPECIFICATIONS	18016
METHOD OF NEBULIZATION:	Ultrasonic Vibration
PARTICLE SIZE:	0.5µm to 5µm
UNIT DIMENSIONS:	3.4" (W) x 5.6" (H) x 2.0" (D)
WEIGHT:	0.51 lbs. or 0.77 lbs. w/battery
MAXIMUM PRESSURE:	2.5MHZ Oscillation Rate
OPERATING PRESSURE:	N/A
LITER FLOW:	N/A
ELECTRICAL REQUIREMENTS:	120 VAC 60Hz
WARRANTY:	3 Years

POPULAR NEB ACCESSORIES AVAILABLE	
18031	Neb Bag, 1/ea
NEBKIT 500	Disposable Neb Kit, 50/cs
18254	Reusable Neb Kit, 10/cs
18090F*	Neb Filter, 10 per bag
*For use with: 18002, 18040-Y, 18040-P, 18070, 18071, 18080, 18090-PG, 18090-PB, 18091-PG and 18091-PB	

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

1P = 1-Pack Available